# VIVASIGHT

By ETVIEW

Instructions for Use

## **VivaSight-SL**

For use by trained personnel only

---

**For use only with Ambu® aView™ monitor**

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## VIVASIGHT-SL

**STERILE (ETO) unless package has been opened or damaged.**

**DO NOT REUSE.**

Federal (US) law restricts this device to sale by or on the order of a physician.

**DESCRIPTION**

*Note: VivaSight-SL is the brand name for the ETVIEW Tracheoscopic Ventilation Tube (TVT™).*

VivaSight-SL is a single-use, cuffed PVC endotracheal tube with an embedded video imaging device and light source at its tip, and an integrated single use video/power cable (adapter cable) with connector. VivaSight-SL displays images of the airway onto the Ambu® aView™ monitor for as long as the device remains in place during intubation.

VivaSight-SL has two Murphy eyes flanking the 45° bevel tip. In addition to the cuffed inflation tube, VivaSight-SL has an injection port for imaging lens cleaning that leads to two lumens running along the tube's wall and opens distally to the imaging lens. The port is used for cleaning the lens with air, saline, or other appropriate solutions.

**INDICATIONS FOR USE**

The ETVIEW Tracheoscopic Ventilation Tube (TVT™) is intended for intubation procedures. The ETVIEW Tracheoscopic Ventilation Tube (TVT™) is indicated for use as a temporary artificial airway in adults requiring mechanical ventilation. It is intended for oral intubations.

The TVT™ System is indicated for viewing during non-difficult and difficult intubation procedures, for verifying endotracheal tube and endobronchial blocker placement and repositioning, for viewing during suctioning and for general inspection of the airway.

**CONTRAINDICATIONS**

Do not use laser equipment in the immediate proximity of the VivaSight-SL.

**ACCESSORIES**

Stylet.

**WARNINGS**

- VivaSight-SL is a single use product. Do not reuse: After use, handle and dispose according to local regulations for disposal of infected medical devices. Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.
- Do not attempt to resterilize any part of VivaSight-SL.
- Do not connect VivaSight-SL into any device or outlet other than those supplied by ETVIEW or Ambu.
- The intended performance of the VivaSight-SL is reached only with the disposable devices and accessories provided or specified by Ambu or ETVIEW.
- VivaSight-SL is available in tube sizes of 7 mm, 7.5 mm and 8 mm (internal diameter) and should be used for patients who require one of these sizes.
- Do not use the injection port for suctioning the patient.
- The cuff pressure should not exceed 25 cm H<sub>2</sub>O. Diffusion of a nitrous oxide mixture, oxygen, or air may either increase or decrease the cuff pressure.
- Do not use the camera in the event that inspired gas temperature exceeds 34 °C for an artificially ventilated patient.
- Do not overinflate the cuff. Overinflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or cuff distortion, which may lead to airway blockage.
- Do not use the tube if the cuff is damaged. Care must be taken to avoid

damaging the cuff during intubation.

- Do not use the tube if the cuff is damaged. Various bony anatomical structures such as teeth, or any intubation aid with sharp surfaces in the intubation path can damage cuff integrity. Care must be taken to avoid damaging the thin-walled cuff during intubation, which may result in the patient requiring traumatic extubation and reintubation.
- The cuff needs to be completely deflated prior to repositioning the tube. Movement of the tube with an inflated cuff may result in damage to the cuff or in patient injury, which may result in medical interventions.
- If the tube is lubricated prior to intubation, it is essential to verify that the lubricant has not entered or occluded the tube lumen or cuff-filling system, thereby preventing ventilation or damaging the cuff. Ventilation and proper functioning of the cuff may be impaired.
- Reliance on the graduated black depth marks on the tube should not substitute for expert judgment. The user should be aware of anatomical variations, including the length of the airway. Intubation and extubation should be performed following currently accepted medical techniques.
- VivaSight-SL's location should be verified every time the patient is moved. Should extreme flexion of the head (chin-to-chest) or movement of the patient (e.g., to lateral or prone positions) occur after intubation, ensure that VivaSight-SL remains in place.
- Do not use an intubation stylet other than provided with the VivaSight-SL.
- Electronic equipment and the VivaSight-SL system may affect the normal function of each other. If the VivaSight-SL system is used adjacent to or stacked with other equipment, observe and verify normal operation of both the VivaSight-SL system and the other electronic equipment prior to using it. It may be necessary to

adopt procedures for mitigation, such as reorientation or relocation of the equipment or shielding of the room in which it is used. Consult the tables in appendix 1 (English version) for guidance in placing the VivaSight-SL system.

- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- The VivaSight-SL is not to be used when delivering flammable anaesthetics to the patient. This could potentially cause patient injury.
- Be careful to check whether the image on the screen is a live image or a recorded image.
- The surface temperature of the VivaSight-SL is likely to reach above 43 °C when the tube is placed outside the patient. Therefore switch off the monitor after the preuse test and switch it back on immediately before use, to prevent the risk of having an impact on the tissue. The surface temperature on the Adapter Cable is likely to reach between 41 °C and 43 °C.
- Prior to use, check that the pouch seal is intact and there are no impurities or damage on the product such as rough surfaces, sharp edges or protrusions which may harm the patient.

## CAUTIONS

- Use of HF equipment in the immediate proximity of the VivaSight SL may lead to impaired image quality.
- To ensure secure placement once the connector has been loosened or removed, the connector and the corresponding part of the tube shall be cleaned with ethanol prior to reconnecting the connector.

- Adapter Cable and connectors shall be kept away from patient face to avoid facial harm.
- VivaSight-SL is intended for use by trained personnel only.
- Use the same size VivaSight-SL as the required standard tube size. Expert clinical judgment should be used in choosing the suitable tracheal tube size for each patient.
- Before use, always check for compatibility between the VivaSight-SL and external instruments (e.g. endobronchial blocker, bronchoscope, suction catheter).
- Do not use VivaSight-SL, suction catheter and endobronchial blocker at the same time. Remove the endobronchial blocker from VivaSight-SL and then insert the suction catheter.
- Do not immerse the tube in liquid.
- The use of lidocaine topical aerosols has been associated with the formation of pinholes in PVC cuffs\*. To prevent cuff leaks, expert clinical judgment is required when using lidocaine topical aerosols.
- Use only water-soluble lubricants. Other lubricants might affect the cuffs.
- Avoid applying lubricants to the front of the camera or the inside of the tube, due to risk of affecting ventilation.
- Use a cuff pressure gauge to help monitor and adjust the cuff pressure. Inflation of the cuff by “feel” alone or by using a measured volume of air is not recommended as resistance is an unreliable guide during inflation.
- Syringes, stopcocks, or other devices should not be left on the inflation line connector.
- Do not connect VivaSight-SL to any external monitor equipment other than aView.

\* Jayasuriya KD, Watson WF: “P.V.C. cuffs and Lidocaine-based aerosol”; *Br J Anaesth.* 1981 Dec; 53 (12): 1368.

## COMPATIBILITY

When VivaSight-SL is used with external instruments (e.g. endobronchial blocker, bronchoscope, suction catheter), use the following table of VivaSight-SL effective inside diameter (defined by ISO 16628) measurements as guidance for choosing the appropriate diameter of the external instrument that will be inserted through VivaSight-SL.

Product	Effective inside diameter
VivaSight-SL 7.0 mm	4.2 mm
VivaSight-SL 7.5 mm	5.2 mm
VivaSight-SL 8.0 mm	5.4 mm

## MRI CONDITIONS

Non-clinical testing has demonstrated that VivaSight-SL is MRI conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 tesla or less.
- Spatial gradient field of 720 gauss/cm or less.
- Maximum specific absorption rate (SAR) of 4 W/kg for 15 minutes.

Under these conditions, VivaSight-SL produces a maximum temperature rise of 0.6 °C.

During MRI scanning, disconnect VivaSight-SL from the single use adapter cable and aView™ monitor, as it may be subject to movement caused by the magnetic field.

MRI image quality may be compromised if the area of interest is in the exact area or relatively close to the position of the camera, VivaSight-SL connector, or cuff inflation port.

## STORAGE AND TRANSPORTATION

- Store and transport VivaSight-SL at temperatures between 0 °C and 42 °C, relative humidity between 10-100% and atmospheric pressure between 80-109 kPa.

- Store in a dry, cool and dark place.

## ADVERSE EVENTS

Adverse events associated with the use of VivaSight-SL are the same as those of standard endotracheal tubes. The most common are laryngospasm; vocal cord paralysis; injuries to lips, gums, tongue, teeth; and aspiration of gastric contents. Consult scientific literature for specific adverse reaction information.

Adverse events associated with use of standard endotracheal tubes, are fracture or dislocation of the cervical spine, endobronchial or esophageal intubation, perforation of the trachea or esophagus, and failed intubation.

## TECHNICAL SPECIFICATIONS

The VivaSight-SL is an electrical safety defibrillation-proof type BF applied part.

- Imaging sensor: CMOS
- Resolution: CIF 320 x 240
- Video format: NTSC Composite Video Baseband Signal (CVBS)
- Light source: 2 LEDs (integrated)
- FOV: Horizontal ~76°, Vertical ~56°, Diagonal ~100°
- Focusing range: 12 mm–60 mm
- Power: up to 36 mA @ DC 5V
- Operating ambient temperature: 10-37 °C (50-98 °F)
- Operating relative humidity: 30-75%
- Operating atmospheric pressure: 80-109 kPa
- Operating altitude: ≤2000 m

## INTUBATION RECOMMENDATIONS

Intubation with VivaSight-SL is performed according to currently accepted medical techniques.

The tube position should be inside the trachea above the main carina.

- Laryngoscope: With VivaSight-SL, use of the laryngoscope may be gentler. It can be used to just hold the tongue.

- Stylet: The accompanying stylet is used to make VivaSight-SL more rigid, for use with or without laryngoscope.
- Oral: Approved for oral intubation procedures.

## PRE-INTUBATION PREPARATION

### CHOOSING VIVASIGHT-SL SIZE

- Use the same size VivaSight-SL as standard required endotracheal tube size.
- Use expert clinical judgment when choosing tracheal tube size for each patient.

## CONNECTING THE TUBE, VERIFICATION STEPS & PREPARATION FOR INTUBATION

1. Connect VivaSight-SL to the single use adapter cable. Be careful to align the white arrow on the single use adapter cable with the arrow on the VivaSight-SL tube round connector.
2. Connect the single use adapter cable to the aView™ monitor, be careful to align the arrow on the cable with the arrow indicator on the aView™ monitor. Verify that all electrical connections are stable and secure. Refer to aView™ monitor instructions for use.
3. Press the power button on top of the aView™ monitor for at least one **second** until a live image appears. Refer to aView™ monitor instructions for use. The aView™ monitor user interface will be ready after approximately one minute, when the hourglass symbol disappears. During this time, proceed with the verification steps.
4. Check that the two LEDs at the tip of the tube are illuminated. If they are not illuminated, replace the tube with another.
5. Verify that there is a stable image on the aView™ monitor (projected from the tip of VivaSight-SL). Then turn off the monitor again until immediately prior to intubation.

6. Test the cuff for integrity by inflating and deflating it completely.
7. Become familiar with the feel of the tube and clear pilot balloon.
8. If applicable, verify that external instruments of appropriate size can be passed through the VivaSight-SL without resistance. There is no guarantee that instruments selected solely using the effective inside diameter will be compatible in combination with the VivaSight-SL.
9. Induce anesthesia, as per protocol. Recommendation: Use a ventilator circuit holder to prevent VivaSight-SL from kinking.
10. Suction the oral cavity and airway area before inserting VivaSight-SL.
11. Prepare the patient, and lubricate the tube (if needed), avoiding the area around the camera lens at the distal end of VivaSight-SL.
12. Intubate the patient and remove the stylet. In case of obstructed or obscured view, see "Cleaning the Lens".

## DEPTH SCALING

*Note: Suggested depth is 2 cm to 3 cm above the main carina. VivaSight-SL has been used at greater depths for one lung ventilation, but keep in mind that VivaSight-SL's electronic wire is outside the protective lumen at 28 cm to 30 cm.*

Like with most intubation tubes, a series of depth marks appear on the side of VivaSight-SL.

Depth marks range from 12 cm to 28 cm, in 2 cm increments.

To determine a specific depth, gauge from the carina:

1. Insert VivaSight-SL's camera up to the depth of the carina.
2. Look at VivaSight-SL's depth marks on the side of the tube.
3. Pull VivaSight-SL back from the carina to the desired depth.

## VISIBILITY DURING PROCEDURE

In the unlikely event that the image is not shown during a procedure, disconnect the VivaSight-SL from the single use adapter cable to the aView™ monitor and continue to use VivaSight-SL as a standard airway tube. **Do not attempt to fix the connection or replace the aView™ monitor.** Turn off the aView™ monitor and dispose of the single use adapter cable.

## CLEANING THE LENS

If VivaSight-SL's lens becomes soiled or obscured by secretions, it can be cleaned by injecting air, prescribed liquid medications, or saline, where permitted by institutional policy, into the injection port of the rinsing tube.

Recommended lens cleaning method:

1. Inject 20 cc air into the injection tube port (RED), and then check image clarity.
2. If the image is not sufficiently clear, connect a 10 cc syringe filled with 2 cc of Saline, to the flushing port and push the plunger. Perform this step twice (total of 4 cc of Saline).
3. Fill a 10 cc syringe with 5-10 cc of air, connect it to the flushing port and push the plunger. Perform this step twice.
4. Connect an empty 10 cc syringe, to the flushing port and suck the Saline. Perform this step twice.

## ADDITIONAL NOTES

- Continuous viewing of the airway enables real-time verification of the tube position, early detection of adverse airway events, and secretion management under visual guidance.
- Validation: During static situations where VivaSight-SL is in place, and the patient is not being moved, validation that the system is operating must be performed periodically by observing breathing movements and cardiac pumping movement in the image on the aView™ monitor.

- The device usage time function on the aView™ monitor is not available for VivaSight-SL. For this device, the device usage time will always display "00:00:00". This will not affect the functioning of the system.

## EXTUBATION

- Deflate the cuff completely, as per protocol.
- Extubation is performed slowly, while examining the airway on the aView™ monitor one last time (document or record, as necessary).

## DISCONNECTING VIVASIGHT-SL

1. Disconnect the single use adapter cable from the aView™ monitor.
2. Turn off the aView™ monitor by pressing the power button for at least two seconds.
3. Do not use the same VivaSight-SL for more than 14 days and do not turn the VivaSight-SL camera ON for more than 8 hours in total. During long term use of the VivaSight-SL, use the depth marks on the VivaSight-SL to monitor the tube position.

## DISPOSAL


Discard the single use VivaSight-SL tube, stylet and adapter cable in a suitable biohazard receptacle in accordance with local regulations.

For disposal instructions for the aView™ monitor, refer to aView™ monitor instructions for use.

## MANUFACTURED BY:

ETView Ltd.  
 Catom 2 Street  
 Misgav Business Park  
 M.P. Misgav 2017900, Israel  
 www.etview.com  
**LB10016 REV03**

 0483	
 Manufacture Date	 Expiry Date
	Consult instructions for use
	Do not use if package is damaged
	Do not reuse
<b>STERILE EO</b>	Sterilized Using Ethylene Oxide
	Non-sterile (single use adapter cable only)
<b>Rx ONLY</b>	US Federal law restricts this device to sale by or on the order of a physician
	Defibrillation-proof type BF applied part
 MR Conditional	 Caution
	This product is not made with phthalates
 42°C / 107°F 0°C / 32°F	Temperature for transportation and storage
 100% 10%	Humidity for transportation and storage
 102kPa 80kPa	Pressure for transportation and storage
	 Manufacturer
<b>EC REP</b>	Authorized Representative in the European Community
	Waste Bin Symbol
	This product is not made with natural rubber latex
<b>REF</b>	Reference Number
<b>LOT</b>	Lot Number, Batch Code

## APPENDIX 1. ELECTROMAGNETIC COMPATIBILITY

Like other electrical medical equipment the system requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC) the system must be installed and operated according to the EMC information provided in this manual.

The system has been designed and tested to comply with IEC 60601-1-2 requirements for EMC with other devices.

### Guidance and manufacturer's declaration – electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment Guidance
RF emission CISPR 11	Group 1	The system use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emission IEC/EN 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

## Guidance and manufacturer's declaration – electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 test level	Compliance Level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2, 4, 8, 15 kV air	+/- 8 kV contact +/- 2, 4, 8, 15 kV air	If floors are covered with synthetic material the relative humidity shall be least 30%.
Electrical fast transient / burst IEC 61000-4-4	+/- 2 kV for mains supply lines +/- 1 kV for input / output lines	+/- 2 kV mains supply lines N/A	Mains power quality shall be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 0.5, 1 kV line(s) to line(s) +/- 0.5, 1, 2 kV line(s) to earth	+/- 0.5, 1 kV line(s) to line(s), mains supply +/- 0.5, 1, 2 kV line(s) to earth, mains supply	Mains power quality shall be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% Ut (100% dip in Ut) for 0.5 cycle at 8 angles	100% reduction 0.5 period at 8 angles	Mains power quality shall be that of a typical commercial or hospital environment.  If the use of the system requires continued operation during power mains interruptions the system can be powered by the built in rechargeable battery.
	0% Ut (100% dip in Ut) for 1 cycle	100% reduction for 1 period	
	70% Ut (30% dip in Ut) for 25 cycles	30% reduction for 25 periods	
	0% Ut (100% dip in Ut) for 5 sec.	100% reduction for 5 sec.	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: Ut is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration – electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted Radio Frequency IEC 61000-4-6	3 V RMS 0,15 MHz – 80 MHz  6 V RMS in ISM bands  80% AM at 1 kHz	3 V RMS 0,15 MHz – 80 MHz  6 V RMS in ISM bands  80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the system, including its cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated Radio Frequency IEC 61000-4-3	3 V/m 80-2700 MHz  80% AM at 1 kHz  Additional proximity field test frequencies and levels according to 60601-1-2 table 9	3 V/m 80-2700 MHz  80% AM at 1 kHz  Additional proximity field test frequencies and levels according to 60601-1-2 table 9	Recommended separation distance  $d = 1.17 \cdot \sqrt{P}$  $d = 1.17 \cdot \sqrt{P}$ 80 MHz to 800 MHz  $d = 2.33 \cdot \sqrt{P}$ 800 MHz to 2.7 GHz  Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range b).  Interference may occur in the vicinity of equipment marked with the following symbol.  

NOTE 1 At 80 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey shall be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system shall be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.

b) Over the frequency range 150 kHz to 80 MHz, field strengths shall be less than 3 V/m.

### Recommended Separation Distances Between Portable and Mobile RF Communication Equipment and system

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters and the system as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.33\sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.37 m	0.37 m	0.74 m
1	1.17 m	1.17 m	2.33 m
10	3.70 m	3.70 m	7.37 m
100	11.7 m	11.7 m	23.3 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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