

Ambu[®] ResKey



Key Benefits

- Transparent face shield with filter, allowing rescuer to visually check for patient's lip colour and vomitus
- Prevents contact with patient's saliva or blood
- Suitable for both adults and children
- Elastic ear bands fit over patient's ears to keep the mask in place
- Can be carried on a key chain, in a pocket, a small handbag and a first-aid kit
- Made from 100% latexfree material
- The Ambu ResKey has a filtration efficiency of FFP3

Ambu[®] ResKey

The new Ambu ResKey is a proven, easy to use, and inexpensive barrier that prevents direct contact with a victim's mouth, nose and face while performing rescue ventilation.

The Ambu ResKey helps to protect rescuer and increases bystander's motivation for CPR in emergency situations.

Standard resuscitation techniques can be used without modification.

Materials

Product	Material
Mask	Transparent LDPE film
Softcase	EVA and PU
Foil	Transparent LDPE film
Rubber bands	Nylon (36%) and rubber (64%)
Filter	Filtrete™ HEPA Filter

Size

Type	
Unfolded	28.5 x 24 cm
Softcase	4.0 x 1.5 x 0.25 cm
Weight app	15 g



Ambu ResKey in use